

# Coupling of Oligosaccharides to Thermo Scientific Nunc CovaLink NH Modules

This Tech Note describes the use of oligosaccharides in monostructured solid phase recognition studies where the oligosaccharide is covalently immobilized on Thermo Scientific Nunc CovaLink NH Modules.

The interest in carbohydrates has increased. One of the main reasons for this is due to the results obtained through research in cell communication, e.g. between leucocytes and in bacterial adhesion.

From that point of view carbohydrates are remarkable information molecules compared to peptides and nucleotides.

If, for instance, two identical monosaccharides are combined, 11 different disaccharides can be built. Whereas, if two identical


aminoacids are combined, only one dipeptide can be formed, and with four different monosaccharides 35,560 unique tetrasaccharides can be created. It is no wonder then that carbohydrates by nature have been assigned a role in communication within living creatures<sup>1,2</sup>.

The diversity of polysaccharides has, however, been a problem in the search for biologically relevant carbohydrate structures for use in the struggle against infection and cancer. Not only is it difficult to collect enough material for examination, it has also been difficult to perform carbohydrate chemistry.

However, synthetically made oligosaccharides have become a potential tool in the search for active structures. The technology

now has evolved to a level comparable to that employed in peptide and nucleotide synthesis, hence a wide range of synthetic saccharides are now commercially available.

## Synthesis of oligosaccharide utilizing the 2-(Tri-methylsilyl)-ethyl (TMSEt) group for anomeric protection.




R = acryl group or »protected« sugar residues


A = acetyl group

- 1) Typically 3 equivalents of acetic anhydride and 1.1 equivalents of boron trifluoride-etherate in chloroform or toluene at room temperature for 1 hour. Purification by column chromatography<sup>3</sup>.
- 2) Mercaptopropionic acid (2.5-3 equivalents) and boron trifluoride-etherate (1.2 equivalents) in methylene chloride at 0°C<sup>4</sup>. Purification by column chromatography.
- 3) Deacetylation using sodium methoxide in methanol. Neutralization and chromatographic purification normally produces the deblocked sugar in quantitative yield.

## Coupling to CovaLink NH



## Binding of activated lactose on CovaLink NH



**Coupling of sugar to Nunc™ CovaLink™ NH and detection of binding exemplified by the Erythrina corallodendron – lactose interaction** <sup>5</sup>. Lactose activated (A-Lac) according to the method described above was used:

- 1) A 1:2 dilution of A-Lac was made in CovaLink NH wells. 100 µL 7.8 mm A-Lac in H<sub>2</sub>O were added to wells in column 1. 50 µL H<sub>2</sub>O were added to the rest of the wells, 50 µL were transferred from column 1 to column 2, etc. After mixing, 50 µL were disposed of from column 12. Finally NHS/EDC in aqueous solution (25 µg each/well) was added to all wells (50 µL/well).
- 2) The plate was incubated on shaker at RT for 1.5 hour.
- 3) Blocked using PBS + 0.5% BSA overnight at 4°C (200 µL/well).
- 4) Rinsed with PBS + 0.5% BSA + 0.05% Triton X-100, 3 times 200 µL/well.
- 5) Biotinylated *Erythrina corallodendron* (10 µg/mL PBS + 0.5% BSA) corresponding to 0.36 µm regarding binding sites (100 µL/well).
- 6) Incubated for 2 hours at RT.
- 7) Rinsed using CovaBuffer (PBS + 2 M NaCl + 65 mm MgSO<sub>4</sub>·2 H<sub>2</sub>O + 0.05% Tween 20) 3 times 200 µL/well.
- 8) Horseradish peroxidase conjugated Avidine added (1:1000 v/v in PBS + 0.5% BSA) (100 µL/well).
- 9) Incubated for 1 hour at RT.
- 10) Rinsed using CovaBuffer 3 times 200 µL/well.
- 11) Substrate (H<sub>2</sub>O<sub>2</sub>, OPD; 100 µL/well) 5 min. at RT; the reaction stopped by adding 1M H<sub>2</sub>SO<sub>4</sub> (100 µL/well).
- 12) The results read at 492 nm.

## References

- 1: Nathan Sharon and Halina Lis. Carbohydrates in Cell Recognition. Scientific American, January 1993, 74-81.
- 2: Nathan Sharon and Halina Lis. Lectins as Cell Recognition Molecules. Science, 246, (1989), 227-234.
- 3: Karl Jansson, Stefan Ahlfors, Torbjörn Frejd, Jan Kihlberg, Göran Magnusson, Jan Dahmén, Ghazi Noorl and Kristina Stenvall. 2-(Trimethylsilyl)ethyl Glycosides. 3. Synthesis, Anomeric Deblocking, and Transformation into 1,2-trans 1-0-acyl Sugars. J. Org. Chem., 53, (1988), 5629-5647.
- 4: Mikael Elofsson, Björn Walse and Jan Kihlberg. Building Blocks for Glycopeptide Synthesis: Glycosylation of 3-Mercaptopropionic Acid and Fmoc Amino Acids with Unprotected Carboxyl Groups. Tetrahedron Lett., 32, (1991) 7613-7616.
- 5: N. Gilboa-Garber and L. Mizrahi. A New Mitogenic D-galactosephilic Lectin Isolated from Seeds of the Coral-tree *Erythrina Corallodendron*. Comparison with Glycine Max (Soybean) and *Pseudomonas Aeruginosa* Lectins. Can. J. Biochem., 59, (1981) 315-320.

## Acknowledgement:

This Tech Note is based on work done by Tomas Ekberg, M.Sc., Organic Chemistry 2, Lund Institute of Technology, University of Lund, Sweden.

[www.thermoscientific.com/oemdiagnostics](http://www.thermoscientific.com/oemdiagnostics)

© 2010 Thermo Fisher Scientific Inc. All rights reserved. "Triton" is a registered trademark of Dow Chemical Company, and "Tween" is a registered trademark of Uniqema Americas. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

## Austria

+43 1 801 40 0

## Belgium

+32 53 73 42 41

## China

+86 21 68654588

## Denmark

+45 4631 2000

## France

+33 2 2803 2180

## Germany

+49 6184 90 6940

## India

+91 22 6716 2200

## Italy

+39 02 02 95059 or  
434-254-375

## Japan

+81 3 3816 3355

## Netherlands

+31 78 571 4440

## Nordic/Baltic countries

+358 9 329 100

## North America

+1 585-586-8800

## Russia/CIS

+7 (812) 703 42 15

## Spain/Portugal

+34 93 223 09 18

## South America

+1 585 899 7298

## Switzerland

+41 44 454 12 12

## UK/Ireland

+44 870 609 9203

## Other Asian countries

+852 2885 4613

## Countries not listed

+49 6184 90 6940 or  
+33 2 2803 2180

TILSPNUNCTN09 0410