


# C8 Black Thermo Scientific Nunc LockWell FluoroNunc MaxiSorp and PolySorp for Fluorescence Detection

The new black Thermo Scientific Nunc LockWell format with MaxiSorp and PolySorp wells have a recommended maximum volume of 350  $\mu$ L. These breakable strips have letters and notches on each well for easy identification of individual wells. The strips are suitable for all commonly used automated equipment. The dense pigmentation of black Nunc LockWell modules minimizes background fluorescence. MaxiSorp™ is optimized for binding IgG (antibodies) and PolySorp™ for binding more hydrophobic molecules.

In this study we compare the performance of the black LockWell modules with similar products from other large suppliers: Competitors A, B and B low volume (total volume of 205  $\mu$ L/well). Plate uniformity and binding capability are determined by monitoring the reaction of immobilized horseradish peroxidase (HRP) from coating with a mixture of rabbit IgG and HRP conjugated rabbit anti-mouse IgG. Three of each plate type were tested on three different days.

When comparing high binding surfaces, the highest binding capability and lowest relative standard deviation was found to be LockWell MaxiSorp and Competitor A (using coating volumes of 150  $\mu$ L/well) (Fig. 1). Comparing the LockWell format to Competitor B low volume (using a coating volume of 100  $\mu$ L/well) a very low relative standard deviation and high binding capability was



**Fig. 1.** Fluorescence intensity (antibody binding capability) measured after performing antibody binding assay on different medium binding surfaces. Coating volume of Competitor A, B and MaxiSorp modules was 150  $\mu$ L/well, and coating volume of PolySorp and competitor B low volume was 100  $\mu$ L/well. Mean CV for tested module plates is shown above the respective binding capability column.

found for the LockWell format. The binding capability of the LockWell MaxiSorp surface can easily be increased using a higher coating volume.

## Assay

- Coating overnight at room temperature with antibody mixture (100  $\mu$ L/well using low volume or 150  $\mu$ L/well using standard 96 format).
- Wash three times with washing buffer.
- Addition of HPPA substrate (100  $\mu$ L/well using low volume or 150  $\mu$ L/well using standard 96 format).
- Addition of 50  $\mu$ L/well sodium hydroxide after 14 minutes. The fluorescence was measured on EnVision 2101 using optimized fluorescence protocols with filter sets 340/405 and 485/535 nm.

## Reagents

Antibody mixture consisting of 65 ng/mL HRP conjugated rabbit anti-mouse IgG and 10  $\mu$ g/mL rabbit IgG, diluted in 0.05 M sodium carbonate buffer, pH 9.6.

Washing buffer: 0.15 M PBS, pH 7.2, with 0.05% detergent (Triton X for high binding surfaces and Tween 20 for medium binding surfaces).

Freshly prepared HPPA substrate: 2.5 mm 3-(p-hydroxyphenyl) propionic acid dissolved in 0.1 M TRIS buffer, pH 8.5, and 1.5  $\mu$ L 30% hydrogen peroxide is added to 100 mL substrate.

Binding capability and relative standard deviations between medium binding surfaces was found to be comparable, using coating volumes of 150  $\mu$ L/well (Fig. 2). Comparing the LockWell

format to Competitor B low volume, a very low relative standard deviation is achieved using a coating volume at 100  $\mu\text{L}/\text{well}$ . The binding capability when using the LockWell for low volume coatings is comparable to the Competitor B low volume. The binding capability of the LockWell PolySorp surface can easily be increased using a higher coating volume.

Fluorescence background measurements were performed using excitation and emission filter sets 485/535 and 340/405 nm. The data are not shown, as it was found relatively low for all products, regardless of the surface type and format.

### Conclusion

Data show high binding capability and high uniformity using fluorescence detection on black C8 LockWell FluoroNunc™, demonstrated for different coating volumes by IgG binding assay on both the MaxiSorp and PolySorp surfaces.

## Labo Baza

nowoczesne wyposażenie laboratorium

ul. Topolowa 5

62-002 Jelonek k/Poznanian

tel.: 061 812 57 45


fax: 061 812 57 25

e-mail: [biuro@labobaza.pl](mailto:biuro@labobaza.pl)

[www.labobaza.pl](http://www.labobaza.pl)

[www.thermoscientific.com/oemdiagnostics](http://www.thermoscientific.com/oemdiagnostics)

© 2010 Thermo Fisher Scientific Inc. All rights reserved. "Triton" is a registered trademark of Dow Chemical Company; and "Tween" is a registered trademark of Uniqema Americas. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.



**Fig. 2.**

Fluorescence intensity (antibody binding capability) measured after performing antibody binding assay on different high medium surfaces. Coating volume of Competitor A, B and PolySorp modules was 150  $\mu\text{L}/\text{well}$ , and coating volume of PolySorp and competitor B low volume was 100  $\mu\text{L}/\text{well}$ . Mean CV for tested module plates is shown above each relevant binding capability column.

#### Austria

+43 1 801 40 0

#### Belgium

+32 53 73 42 41

#### China

+86 21 68654588

#### Denmark

+45 4631 2000

#### France

+33 2 2803 2180

#### Germany

+49 6184 90 6940

#### India

+91 22 6716 2200

#### Italy

+39 02 02 95059 or  
434-254-375

#### Japan

+81 3 3816 3355

#### Netherlands

+31 76 571 4440

#### Nordic/Baltic countries

+358 9 329 100

#### North America

+1 585-586-8800

#### Russia/CIS

+7 (812) 703 42 15

#### Spain/Portugal

+34 93 223 09 18

#### South America

+1 585 899 7298

#### Switzerland

+41 44 454 12 12

#### UK/Ireland

+44 870 609 9203

#### Other Asian countries

+852 2885 4613

#### Countries not listed

+49 6184 90 6940 or  
+33 2 2803 2180

TILSPNUNCTN59 0510