SONOPULS

Ultraschall-Homogenisatoren • *Ultrasonic homogenizers*

Zubehör / Accessories

Emulgieren

Homogenisieren

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Suspendieren

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Sonochemie

Emulsifying

Homogenizing

Cell disruption

Suspending

Accelerating reactions

Sonochemistry

BANDELIN

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55 years of experience in ultrasound technology

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- **Amplitude control**
- **Pulsation** 10 - 100 % - storable
- **Timer** continuous or timed operation up to 99 min
- **Automatic storage** of previous settings for pulsation and process duration
- **LED-Display** - monitoring of preselected values for amplitude, pulsation and processing time

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Small unit for lab routine

Ready-to-operate basic equipment for volumes from 2 ml to 50 ml with generator GM 2070, ultrasonic converter UW 2070, stepped standard horn SH 70 G and microtip MS 73, 3 mm diameter. HF-output max. 70 W_{eff}

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Standard unit for lab routine

Ready-to-operate basic equipment for volumes from 20 ml to 900 ml with generator GM 2200, ultrasonic converter UW 2200, booster horn SH 213 G and titanium flat tip TT 13 of 13 mm diameter. HF-output max. 200 W_{eff}

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The **AMPLICHRON®-system** guarantees a constant amplitude independently from changing conditions within the sample. It ensures reproducible results for process validation. Settings within a range of 10 to 100 %. Verification of actual value at the display.

Permanent control of ultrasound irradiation as well as indication of wear of the probe.

Pulsation

permits heat dissipation when processing sensitive samples. The adjustable pulsation allows cooling during rest intervals.

Built-in timer

Process duration storable. Indication of elapsed time in continuous operation or remaining time in countdown mode.

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- **Amplitude control** or alternatively power control
- **Pulsation:** operation intervals from 0,1 s up to 600 s, rest intervals from 0,2 s up to 600 s
- **Timer** up to 9 h: 59 min: 59 s
- **Display of energy (kJ)** delivered to the sample
- **9 user programs**
- **Monitoring of allowed probes**
- **Automatic amplitude limiting** after preselection of probe
- **Indication of temperature** from 0 up to 120 °C, alternatively switching off or warning signal, temperature sensor optionally
- **WINPULS remote control** for process organization with PC over RS-232 and IR interface
- **Lighted LCD display**

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High tech for research

- for small volumes -
Ready-to-operate basic equipment for volumes from 2 ml to 50 ml with generator GM 3100, ultrasonic converter UW 3100, stepped standard horn SH 70 G and microtip MS 73, 3 mm diameter.
HF-output max. 100 W_{eff}

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High tech for research

- for larger volumes -
Ready-to-operate basic equipment for volumes from 20 ml to 900 ml with generator GM 3200, ultrasonic converter UW 3200, booster horn SH 213 G and titanium flat tip TT 13 of 13 mm diameter.
HF-output max. 200 W_{eff}

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Processing frequency - 20 kHz uniform sonic waves

Constant sound radiation- extremely effective.

Switching ON / OFF - easy to handle

either at the generator or directly at the ultrasonic converter via button or remote control.

Compact construction

Space saving design. Ideal for use in mobile laboratories.

Foil keypad

easy to clean and user-friendly

Fail-safe during continuous operation and idling

RFI-proofed and CE-marked, also as medical device compliant to the directive for in-vitro diagnostics 98/79/EG

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Processing time

up to 99 min indication of remaining time or continuous operation with indication of elapsed time

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Built-in parameter memory

automatic storage of previous settings for duration of ultrasonic process and pulsation

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
Pulsation 10 - 100 %. Safe reproduction through value storage and display control.

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Amplitude indicator measured value
in %

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Adjustable amplitude
linear adjustment of amplitude (μm) to set-point

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Mode key
programming and storage of time and pulsation mode

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Reset key
restart / selection of operating mode

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Display of measured value
amplitude (Am) in % or power (Pw) in W

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
Processing time

up to 9 h 59 min 59 s indication of remaining time or continuous operation with indication of elapsed time

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All-in time
(only in countdown operation)

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How to select the proper unit

Power output in watt is not the sole criterion for selecting an ultrasonic homogenizer. This value only indicates the power of the HF-generator but not the energy delivered to the sample. The amplitude at the radiating surface of the probe is the determining factor for the evaluation of the irradiation result while taking into consideration the volume of the sample.

SONOPULS ultrasonic homogenizers produce with the same electric power higher amplitudes than other models in the market thanks to an ideal integration of all components.

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Ultrasonic homogenizers are used in laboratories, hospitals and industry for scientific experiments, analysis as well as in pilot or small-lot production. Here are some examples showing the vast variety of applications for ultrasonic homogenizers:

Chemistry and sonochemistry

- *acceleration of chemical reactions*
- *increase of reactivity of metallic interfaces*
- *destroying of highly-molecular compounds - depolymerisation*
- *catalytic effect during reaction of powders: increase of the catalytic effect of nickel powder by more than 100.000 times*
- *cracking of alkanes into smaller fragments at room temperature*
- *increasing of reaction rate during Ullmann reaction, grignard reaction, Wurtz reaction*

Biology / medical range


- *fast disruption, e. g. of Acetobacter xylinum, Actinomyces, Bacillus subtilis, Candida albicans, Erythrocytes, Escherichia coli, Eukaryotic cells, H. Ep. 2-cells, Helicobacter pylori, Saccharomyces cerevisiae, Lactobacillus, Staphylococcus aureus, Streptococcus*
- *isolation of heads from spermatozoes*
- *simplification of paternity test of stromafree hemolysate, reducing of processing time from 30 min to 5 s*
- *breaking of DNA-chains*
- *separation of enzymes from mitochondria*
- *acceleration of production of misfolding prions in medicine research*
- *cleaning of radioactive contaminated seeds for non-leakage-proof*
- *producing of homogenous infecting solutions; processing of poxviridae*
- *cell disruption for blot Fleck*

Techniques and analysis

- *production of stable and finest emulsions*
- *preparation of soil samples for determination of mineral-oil like hydrocarbons, polycyclic aromatic hydrocarbons (PAH)*
- *homogenization of waste samples to define possibility of waste disposal and to characterize pollutants*
- *dispersion of agglomerates, e.g. in colour and lacquer production*
- *dispersing of samples for determination of grain size*
- *degassing of liquids*

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Analysis

Dispersing of solid particles

Aim: grain size analysis
Quantity: 50 - 100 ml
Approx. time: 2 - 5 min
Unit: HD 2200/3200 with KE 76

Waste water samples

Aim: Homogenizing for determination of harmful substances, e.g. mineral oil, grease AOX in industrial and butcher's waste water
Quantity: 250 ml
Approx. time: 5 - 10 min
Unit: HD 2200/3200 with TT 13, for taller vessels with VS 70 T

Soil samples

Aim: Extraction for determination of pH value, Mg, K, P, N – contents for recommendation of fertilizer / determination of radio nucleides to control radioactivity in the environment (milk researchment)
Quantity: 50-100 ml / 100-150 ml
Approx. time: a few seconds
Unit: HD 2200/3200 with KE 76 / VS 70 T

Meat and sausages samples

Aim: Homogenizing for determination of nitrates
Quantity: 100 ml
Approx. time: 3 min
Unit: HD 2200/3200 mit KE 76

Cheese

Aim: Homogenizing for determination of nitrates
Quantity: 10 g cheese in 25 ml water
Approx. time: ½ min
Unit: HD 2200/3200 with KE 76

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Biochemie - Biologie - Medizin

Candida albicans

Aim: Disruption
Quantity: 10 ml
Approx. time: 10 min
Unit: HD 2070/3100 with MS 73

Escherichia coli

Aim: Disruption for proteine lay off
Quantity: 10 ml
Time: 5 - 10 min
Unit: HD 2070/3100 with MS 73 or HD 2200 with MS 73

Eucaryotic cells

Aim: Disruption for proteine lay off
Quantity: 1,5 ml
Approx. time: 1 min
Unit: HD 2200/3200 with BR 30 and EH 3, cooling necessary

Heart muscle tissue

Aim: Homogenizing
Quantity: 1,5 ml
Approx. time: 4 min
Unit: HD 2200/3200 with MS 72, cooling necessary

Insect cells

Aim: Disruption for proteine lay off
Quantity: 20 - 50 ml
Approx. time: app. 25 sec, pulsed
Unit: HD 2070/3100 with MS 73 and RZ 2

Kidney tissue

Aim: Homogenizing
Quantity: 1,5 ml
Approx. time: 40 sec
Unit: HD 2200/3200 with MS 72, cooling necessary

Liver tissue

Aim: Homogenizing for molecular genetic tests
Quantity: 1,5 ml
Approx. time: 1½ min
Unit: HD 2200/3200 with MS 72, cooling necessary

Liposomes

Aim: Producing of small unilamellar phospholipid vesicles
Quantity: 20 ml
Approx. time: 10 - 15 min
Unit: HD 2070/3100 with TT 13, cooling necessary

Lymphocytes

Aim: Disruption
Quantity: 50 µl - 2 ml
Approx. time: 1 - 5 min
Unit: HD 2070/3100 with BR 30 and EH 3

Saccharomyces cerevisiae

Aim: Disruption
Quantity: 20 ml
Approx. time: 15 min
Unit: HD 2200/3200 with KE 76, addition of glass beads to accelerate process, cooling necessary

Staphylococcus aureus

Aim: Disruption
Quantity: 10 ml
Approx. time: 10 min
Unit: HD 2070/3100 with MS 73

Streptococcus

Aim: Disruption
Quantity: 10 ml
Approx. time: 8 - 10 min
Unit: HD 2200/3200 with MS 73

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Probes made of titanium alloy (Ti-6Al-4V) transmit mechanical longitudinal waves into the sample. They are thermo-resistant, can be treated in autoclaves and are resistant to corrosive media. Sample volume, diameter of the processing vessel and the required amplitude determine the selection of the unit and the type of probe. The higher the amplitude, the more intense the sonication.

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VS recommended for long-necked vessels.
VS 70 between SH 70 G / 213 G and TT 13
VS 200 between SH 225 G and TT 25

Typ / type **VS 70**
 für / for HD 2070/2200/3100/3200
 Best.-Nr. / code no. 500

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
GS for connection to HD 2070/3100 with **special horn SH 70 GQ**. For application in food analysis, pharmacy or environmental analysis. No intrusion of metal particles and boron compounds - ideal for trace analysis. High chemical and temperature shock resistance, no electric conductivity.

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MR for connection to HD 2070/3100 composed of Multison horn MRH 12 and 12 Multison tips MRS 2 (dia. 2 mm) or MRS 3 (dia. 3 mm). For irradiation of samples in microplates and deep well plates. Simultaneous sonication of 12 samples. MRS made of titanium alloy; individually replaceable.

Typ / type **MR 12-2** **MR 12-3**
 MRS Ø / dia. 2 mm 3 mm
 Best.-Nr. / code no. 3626 3633

Wendelsonotrode • Spiral probe - Gebrauchsmuster / registered patent D 203 07 648

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WS 8, dia. 8 mm, with **special horn SH 200 W**, for connection to HD 2200/3200. For gentle use in slim reaction vessels. Lateral irradiation.

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
Standard and booster horns (Ti-6Al-4V) are furnished with a thread for replaceable probes. With exterior thread to connect various vessels.

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Solid standard horn
DH 13 G with diamond coating on the radiating area; lifetime is thirty times longer than usually.

Typ / type	SH 70 G	SH 213 G	SH 219 G	SH 225 G	DH 13 G
für / for HD	2070	2200	2200	2200	2070/2200
	3100	3200	3200	3200	3100/3200
Best.-Nr. / code no.	486	527	600	543	403

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FZ 5 G und FZ 7 G to prepare stable mixtures of non-mixable or hardly mixable liquids (oil-in-water) by direct intrusion of pre-mixed samples into the cavitation field. In combination with flow-through cell DG 4 G the continuous treatment of 2 different media and parallel tempering is possible. Material: Ti-6Al-4V

Typ / type	FZ 5 G	FZ 7 G
für / for HD	2070/3100	2200/3200
Best.-Nr. / code no.	490	452

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Sleeve adapters made of PTFE for tight mounting on standard ground glass vessels.
NA 29 G for NS 29/32 for SH 70/213 G
NA 45 G for NS 45/40 for SH 70/213/219/225G
GA 3 G threaded adapter made of stainless steel with external thread M 40 x 1 for SH 70/213/219/225 G

Typ / type	NA 29 G	NA 45 G	GA 3 G
für / for HD	2070/2200	2070/2200	2070/2200
	3100/3200	3100/3200	3100/3200
Best.-Nr. / code no.	540	487	473

Titaneller titanium flat tips			Quarzglas-Sonotroden silica glass probes						Wendelsonotrode spiral probe
TT 13	TT 19	TT 25	GS 6	GS 6 L	GS 13	GS 13 L	GS 18	GS 18 L	WS 8
497	491	532	024	048	028	050	040	054	3617
13	19	25	6		13		18		8
5	5	6	145	290	145	290	145	290	215
SH 70 G	-	-	SH 70 GQ		SH 70 GQ		SH 70 GQ		-
SH 213 G	SH 219 G	SH 225 G	-		-		-		SH 200 W
78 / 93	- / -	- / -	12,5 / 13		13 / 13		13 / 13		- / -
149 / 165	73 / 81	48 / 53	- / -		- / -		- / -		12 / 12
10 - 200	-	-	2 - 100		25 - 200		25 - 500		-
20 - 900	25 - 900	30 - 1000	-						8 - 15
17	23	29	10		17		22		10

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Processing vessel, stainless steel

DG 4 G for high-volume flow-through processing like emulsifying, dispersing or homogenizing, up to 30 l/h,

KG 4 G, closed vessel with cooling jacket. Processing volume about 65 ml.

Typ / type	KG 4 G
für / for HD	2070/2200/3100/3200
Best.-Nr. / code no.	3609

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Processing vessels made of glass

Cooling vessel KG for sonication of temperature-sensitive samples.

Flow-through vessel DG with cooling jacket for irradiation of larger volumes.

Rosett cell RZ for homogenous and intense circulation of sample.

Suslick cell SZ with 3 sidearms for introduction of gas or measuring probes.

Typ / type	KG 3	KG 5	DG 3	DG 5
für / for HD	2070/2200	2200	2070/2200	2200
3100/3200	3200	3100/3200	3200	
Volumen / volume	15 ml	80 ml	max. 5,6 l/h	max. 30 l/h
innen / interior diameter	20 mm	35 mm	20 mm	35 mm
Tiefe / height	65 mm	95 mm	65 mm	95 mm
Best.-Nr. / code no.	536	481	538	482

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Cup Horn BB for indirect irradiation with high energy. No contamination from cavitation erosion. No aerosoling with pathogenic or hazardous materials.

BB 6: temperature regulation with circulating liquid coolants.

Microtube holder EH 6 for use in BB 6. Up to 6 samples can be treated simultaneously. Pressure plate holds tubes in place. No floating of cups.

Cup booster BR 30 for high-intensive irradiation of smallest and sensitive sample volumes, e. g. cleaning of radioactively contaminated seeds, cell disruption or treatment of pathogens and toxic material. No cross-contamination or sample loss. Flow-through cooling liquid avoids rise of temperature.

Inset basket BK 30 for intensive cleaning of small parts in BR 30.

Microtube holder EH 3 for use with BR 30. Up to 3 samples can be treated simultaneously. 2 exchangeable discs with diameters 8,5 or 11,5 mm.

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Stainless steel stand HG 5
with lab clamp, for safe fixing of ultrasonic converter.

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Sound proof box
reduces noise level. Inlets at rear side for gas supply or flow-through irradiation. Clear acrylic glass door permits observation of process. LS 6 with height adjustable stand. **New:** For indirect processing with BR 30, BB 6 or BB 2 G turn the box by 180°. Ultrasonic converter UW is fixed safely by special clamp. **Support UG 6** is required.

Typ / type für / for HD	LS 4 2070/2200 3100/3200	LS 6 2070/2200 3100/3200	UG 6 2070/2200 3100/3200
Dämpfung / noise reduction in dB-AU	10	20	
Best.-Nr. / code no.	416	3615	3616

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Distance tube AH 6
to be directly clamped into the closure of the LS 6, for direct sonication with long probes (MS 72/73, KE 76, VS 70 with TT 13, VS 200 with TT 25, VS 70 T, VS 200 T and GS..)

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
Foot switch remote control
For easy switching ON/OFF of the HF-generator. With 3 m cable.

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
Temperature sensor
for measuring of sample temperature from 0 up to 120 °C.

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For process control with PC for operation systems WIN 2000 and WIN XP with different additional functions like test logging and comfortable data storage. (up to 99 storages). Set composed of WINPULS®-Software and infrared adapter IR 1 for interface RS 232.

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Best.-Nr. / code no. 8067

Applications

- *biotechnology and sonochemistry*
- *homogenizing of waste water for environmental analysis*
- *emulsifying in cosmetic industry, photolithography etc.*
- *suspending*

Features

- *efficient omnidirectional ultrasonic radiation*
- *flow-through rate max. 50 l/min*
- *high power density 350 W/l*
- *ultrasonic frequency 40 kHz*
- *microprocessor-controlled*
- *power output constancy $\pm 2\%$*
- *modular extension up to 2000 W*

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Further documentation on accessories and other SONOREACTORS on request.

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- *standard baths from 1 - 200 litres*
- *compact construction*
- *easy to handle*
- *optionally with heating*
- *integrated timer*
- *from 3 ltr. upwards with drain*
- *cost-saving - reduced cleaning time*
- *intensive, gentle, thorough*
- *large range of accessories*

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BANDELIN *electronic* - specialized in the manufacturing of ultrasonic units, maintains a Quality Management System complying with the requirements of EN ISO 9001/12.2000 and EN ISO 13485/11.2000

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