

**NALGENE®**

**nunc™**

# **RNase/DNase-Free Products and Capabilities**

## **Labo Baza**

**nowoczesne wyposażenie laboratorium**

ul. Topolowa 5

62-002 Jelonek k/Poznania

tel.: 061 812 57 45

fax: 061 812 57 25

e-mail: [biuro@labobaza.pl](mailto:biuro@labobaza.pl)

**[www.labobaza.pl](http://www.labobaza.pl)**

- **Cryobank and Bank-It™ Storage Systems**
- **Cryogenic Vials**
- **Storage Vials**
- **Micro Packaging Vials**
- **Centrifuge Tubes**
- **MicroWell® Plates**
- **Amplification Plates and Tubes**
- **Sealing Tapes**
- **Reagent Reservoirs**

**Nalge Nunc  
International**

## NUNC™ Sealing Tapes for Plates

- Prevents contamination and moisture loss
- Polyolefin tape has good optical clarity and low autofluorescence
- Aluminum tapes are pierceable

NUNC™ Cat. No.	Material	Color	Note	Units per pack/case
232702	Polyolefin	Clear	General use, PCR	100/100
276014	Aluminum	Silver	General use, storage	100/100
232698	Aluminum	Silver	PCR, Cold Storage, 96 well plates	100/100
232699	Aluminum	Silver	PCR, Cold Storage, 384 well plates	100/100



## NALGENE® and NUNC™ Sterile Disposable Polypropylene Robotic Reservoirs


- 300ml volume in a MicroWell® plate footprint; 86 x 128 x 43 mm
- Convoluted bottom design minimizes dead volume and has baffles to reduce splashing
- Automation-compatible

NALGENE® Cat. No.	NUNC™ Cat. No.	Bottom	Units per pack/case
1200-1301	370906	Flat	4/40
1200-2301	370908	Convoluted	4/40

## NUNC™ 15 ml Polypropylene Centrifuge Tubes with HDPE Closure

- Available in polypropylene for chemical resistance
- Printed volume graduations along length of tube, molded graduations on conical section
- Meets USP Class VI
- Flat cap and frosted area on tube provide writing surfaces for sample identification
- Tube and leakproof cap suitable for freezing
- PP tube withstands temperatures up to 125°C

NUNC™ Cat. No.	Max. RCF tested (xg):	Dimensions, mm	Packaging	Units per pack/case
366036	7800	120 x 17	Racked	50/500
366060	7800	120 x 17	Bulk	50/500



## NUNC™ PPCO Vials with Polyethylene Screw Cap, High Purity Polypropylene Copolymer

- Excellent for centrifugation and storage of PCR reagents, restriction enzymes and other biochemical reagents
- Frosted writing area
- Leakproof closure
- Suitable for storage down to -95°C
- Racks and color-coded cap inserts available

NUNC™ Cat. No.	Max. RCF tested (xg):	Dimensions, mm	Working Vol., ml	Total Vol., ml	Units per pack/case
264262	13,000	49 x 13	0.5	0.9	500/500
264300	13,000	49 x 13	2.0	2.2	500/500

