

nunc™

NALGENE®

Storage Racks for MicroWell® Plates



- Store, categorize and transport MicroWell Plates
- Plastic, aluminum or stainless steel
- Accommodates low profile, standard height and DeepWell® plates
- Freezer and incubator compatible
- Colors and adjustable shelves available

Aluminum Plate Storage Racks in Colors

- Adjustable shelves to fit any plate or tube rack
- Color coded top and sides for easy identification
- Locking bar keeps plates secure
- For chest and upright mechanical freezers or incubators
- Lightweight aluminum
- Autoclavable
- Also excellent for cluster tube racks



NUNC Cat. No.	367001	367002	367003	367004
Capacity (plates)	18 Low Profile 10 Standard 6 DeepWell 1ml 5 DeepWell 2ml	18 Low Profile 10 Standard 6 DeepWell 1ml 5 DeepWell 2ml	18 Low Profile 10 Standard 6 DeepWell 1ml 5 DeepWell 2ml	18 Low Profile 10 Standard 6 DeepWell 1ml 5 DeepWell 2ml
Color	Blue	Natural	Yellow	Red
Configuration	5 - 18 Tall	5 - 18 Tall	5 - 18 Tall	5 - 18 Tall
Rack Size, H x W X D, cm	23.5 x 9.5 x 13.5	23.5 x 9.5 x 13.5	23.5 x 9.5 x 13.5	23.5 x 9.5 x 13.5
Height per shelf, cm	Adjustable	Adjustable	Adjustable	Adjustable
Units per pack/case	1	1	1	1

Stainless Steel Plate Storage Racks



- Durable stainless steel construction
- Retainer clip on each compartment holds plates securely in place
- Sturdy handles on top and sides for easy transport
- Label holder for easy identification
- Holds standard plates or deep well plates
- Two convenient horizontal sizes to fit your freezer or refrigerator

NALGENE Cat. No.	5039-0048	5039-0072
Capacity (plates)	80 Low Profile 48 Standard 16 DeepWell	120 Low Profile 72 Standard 24 DeepWell
Color	Natural	Natural
Configuration	4 Tall x 4 Wide	4 Tall x 6 Wide
Rack Size, H x W X D, cm	26.7 x 48.3 x 17.8	26.7 x 63.5 x 17.8
Height per shelf, cm	5.3	5.3
Units per pack/case	1	1


Labo Baza

nowoczesne wyposażenie laboratorium

ul. Topolowa 5
62-002 Jelonek k/Poznania
tel.: 061 812 57 45
fax: 061 812 57 25
e-mail: biuro@labobaza.pl
www.labobaza.pl

NUNC Storage Plates and Accessories

- 96- and 384-well polypropylene MicroWell plates
- Standard and DeepWell® Plates
- Colors available
- Barcoding available
- Sealing tapes, Lids and Capmats to protect samples



See the NUNC Catalog or our website at www.nuncbrand.com for the full line of products.