



### Compact Shaker KS 15 ↻ or ↔

Reliable, robust shaker for daily use

<b>Motion:</b>	<b>Model A: orbital</b> <b>Model B: reciprocating (horizontal)</b>
<b>Shaking platform:</b>	400 x 300 mm
<b>Max. load:</b>	15 kg
<b>Shaking speed:</b>	30 - 420 rpm
<b>Stroke:</b>	17 mm
<b>Timer:</b>	0 - 120 min / continuous

#### Basic Equipment:

Basic device **without shaking plate**, without rack system

**KS 15 A: No. 6171 000**  
**KS 15 B: No. 6173 000**

#### Rack Systems:

- ① Rack system Combifix KS ..... No. 0052 071
- ② Universal tray KS ..... No. 0051 471
- ③ Incubator hood TH 15 ..... No. 6161 000

Accessories: see page 22

#### Loading Capacity:

250 ml Erlenmeyer flasks: 12/15 pieces  
500 ml Erlenmeyer flasks: 6/9 pieces  
250 ml separating funnels: 4 pieces  
Further loading possibilities: see page 23

#### Technical Data:

Speed range:	30 - 420 rpm, steplessly adjustable
Stroke:	17 mm
Max. load:	15 kg
Timer:	0 - 120 min or continuous operation
Electrical supply:	230 V or 115 V, 50/60 Hz (please indicate in case of order)
Enclosure protection:	IP 21
Heat emission:	approx. 10 W
Ambient temperature:	5° C to 50° C
Relative humidity:	~ 85 %
Dimensions (wxdxh):	515 x 495 x 235 mm
Weight:	22 kg



### Compact Shaker KS 15 control ↻ or ↔

For reproducible shaking tasks

<b>Motion:</b>	<b>Model A control: orbital</b> <b>Model B control: reciprocating (horizontal)</b>
<b>Shaking platform:</b>	400 x 300 mm
<b>Max. load:</b>	15 kg
<b>Shaking speed:</b>	30 - 420 rpm
<b>Stroke:</b>	17 mm
<b>Timer:</b>	1 min - 30 h / continuous

#### Basic Equipment:

Basic device **without shaking plate**, without rack system; with digital display (set values/actual values) and programmable control

**KS 15 A control: No. 6170 000**  
**KS 15 B control: No. 6172 000**

#### Rack Systems:

- ① Rack system Combifix KS ..... No. 0052 071
- ② Universal tray KS ..... No. 0051 471
- ③ Incubator hood TH 15 ..... No. 6161 000

Accessories: see page 22

#### Loading Capacity:

100 ml Erlenmeyer flasks: 20/24 pieces  
1000 ml Erlenmeyer flasks: 4/5 pieces  
Further loading possibilities: see page 23

#### Technical Data:

Speed range:	30 - 420 rpm, steplessly adjustable
Stroke:	17 mm
Max. load:	15 kg
Timer:	1 min - 30 h or continuous operation
Electrical supply:	230 V or 115 V, 50/60 Hz (please indicate in case of order)
Enclosure protection:	IP 21
Heat emission:	approx. 18 W
Ambient temperature:	5° C to 50° C
Relative humidity:	~ 85 %
Dimensions (wxdxh):	515 x 495 x 235 mm
Weight:	22 kg

## COMPACT SHAKERS

Accessories/Rack Systems for Compact Shakers



① ① Rack system Combifix KS  
No. 0052 071



② ② Universal tray KS  
No. 0051 471



③ ③ Incubator hood TH 15  
No. 6161 000

## Labo Baza

nowoczesne wyposażenie laboratorium  
ul. Topolowa 5  
62-002 Jelonek k/Poznań  
tel.: 061 812 57 45  
fax: 061 812 57 25  
e-mail: biuro@labobaza.pl  
[www.labobaza.pl](http://www.labobaza.pl)

Detailed description of the rack systems: see page 18  
TH 15: see page 14